



DERMINE



1.1 ML FINE / LIPS / DEEP

DM
LABORATORIES



Enhance your beauty



DERMINE is the new generation of cross-linked Hyaluronic Acid filler. It is a monophasic filler mono-linked hyaluronic acid, which is ideal for subcutaneous applications. Injectable gel 1ml with 2 needles. Hyaluronic Concentration 24mg/ml. Lidocaine 3mg/ml.



Injectable gel 1ml with 2 needles. Hyaluronic Concentration 24mg/ml. Lidocaine 3mg/ml.



- ULTIMATE Thin Wall needle
- The Ultra Thin Wall needle is thinner wall, which means more smooth injection.



Product Model	Facial Area	Injection Sites Shown in Image	Recommended Injection Depth
FINE	Upper facial area	Eyebrow region, lateral temple, fine periorbital lines	Superficial dermis
DEEP	Midface and jawline	Cheeks, jawline, chin, lateral facial contours	Subcutaneous or supraperiosteal
LIPS	Perioral area	Upper lip, lower lip, perioral lines, nasolabial folds adjacent to the lips	Deep dermis / submucosal

PRODUCT MODELS & INDICATIONS



FINE

silver

Description

- Injection Depth: Superficial Dermis
- Indications: Fine facial lines, superficial wrinkles, perioral lines, forehead fine lines

\$12.5



DEEP

gold

Description

- Injection Depth: Subcutaneous or Supraperiosteal
- Indications: Volume creation, facial contours, cheeks, chin, nose bridge, forehead

\$12



LIPS

black


Description

- Injection Depth: Deep Dermis / Submucosal
- Indications: Lip augmentation, deep facial folds, nasolabial folds, marionette lines


\$12

Manufacturing of

DERMINE



Certifications & Quality Standards



ISO & Quality Management Certifications:

DERMINE is manufactured in compliance with internationally recognized quality management standards, ensuring safety, consistency, and reliability for medical use. The manufacturing facility of DM LABORATORIES Co., Ltd. is certified according to the following standards:

- ISO 13485:2016
- Medical Devices – Quality Management System
- Scope: Manufacture of non-active implantable medical devices
- ISO 9001:2015 / KS Q ISO 9001:2015
- Quality Management System
- Scope: Manufacture of medical devices

These certifications confirm that the manufacturer operates under a controlled and audited quality system in accordance with international regulatory requirements

STORAGE & MANUFACTURER INFORMATION

Storage Conditions

- Store between 1°C and 30°C
- Protect from freezing and direct sunlight

Sterilization

- Product: Steam sterilization
- Needle and syringe: Ethylene oxide sterilization

Manufacturer

DM LABORATORIES Co., Ltd
16th Floor, Gangnam Building
1321-1 Seocho-dong, Seocho-gu, Seoul 137-070

DM LABORATORIES **DERMINE**

 [dermine.kr](https://www.instagram.com/dermine.kr)

 www.dermine.kr

